



DATAcentre Base Station Memory Upgrade

Compatible with DATAcentre

- MCC card memory upgrade
- Up to 10 years on-board storage
- Card formatted to work with DATAcentre software
- Up to 32GB available



Each DATAcentre Base Station comes with a 2GB memory as standard. However depending on the number of points being monitored and the time that data is to be stored on the base station this memory can be upgraded.

Planer specifically format the memory to ensure it will work with the base station software; non Planer supplied memory cards may not work correctly with the system and data may not be stored in the correct format.

Technical Specifications				
Memory Capacity	16 Points	32 Points	64 Points	120 Points
2GB (std)	5.1 years	2.7 years	1.4 years	0.7 years
4GB	10.2 years	5.3 years	2.7 years	1.5 years
8GB	20.3 years	10.7 years	5.5 years	3 years
16GB	40.7 years	21.3 years	10.9 years	5.9 years
32GB	81.4 years	42.6 years	21.8 years	11.7 years

Part Numbers	
FTEM02MC101	4GB
FTEM02MC102	8GB
FTEM02MC103	16GB
FTEM02MC104	32GB