



Supplied with 1.5 m 3 core flat PVC ribbon cable

PT100 Surface Temperature Sensor -20 °C to +90 °C

DATAssure™ D3

FTEP3-DS305-HR

For use with the DATAssure™ D3 family of wireless base stations. Compatible with a full range of Intelligent Wireless Transmitters (supplied separately).

Polyimide insulated surface mount temperature sensor designed to provide a practical method for measuring surface temperatures in the range of -20 °C to +90 °C. The thin design of the sensing element minimises the effects of having the sensor attached to the surface being measured resulting in a fast thermal response. The sensor is terminated to 3 wire flexible flat ribbon cable providing easy mount through door seals.

Common uses: Chilled & Hot Surfaces | Hot Plates | Incubators

This product can be supplied with a calibration certificate.

Technical Specification			Sensor Cable:	1500 mm 3 core PVC ribbon cable
Temperature Range:			-20 °C to +90 °C	
Accuracy:				
	-20 ° to 40 °C	40 ° to +90 °C	Cable Operating Temperature:	-20 °C to +90 °C
Standard	- ±0.5 °C	±0.8 °C	Probe Sheathing:	Polymide Film
Calibrated	- ±0.2 °C	±0.3 °C	On Cable Electronics	
Thermistor:	PT100 Class B		Connector IP Rating:	IP67
Weight:	19 g		Connector:	3 pin male Hirose
			Operating Temperature:	-20 °C to +55 °C
			Electronics Sheathing:	Technomelt OM678
			Electronics & Probe IP Rating:	IP68

Compatible FTEP2-WS Wireless Transmitters

Part No.	Description
FTEP2-WS102-0433	Single Input Wireless Transmitter with Hirose Connector (433 Mhz)
FTEP2-WS103-0433	Dual Input Wireless Transmitter with Hirose Connectors (433 Mhz)
FTEP2-WS102-2400	Single Input Wireless Transmitter with Hirose Connector (2.4 Ghz)
FTEP2-WS103-2400	Dual Input Wireless Transmitter with Hirose Connectors (2.4 Ghz)

Ci070/V1