



<http://www.steris-ast.com>

Certificate of Irradiation

Date Issued: 24-Jun-2021

FR01S12579241-1-1

This is to certify that Synergy Health Marseille, a STERIS Company has where appropriate delivered an irradiation process in accordance with the current certified standards:

EN ISO 11137-1 Sterilisation of Health Care Products
EN ISO 13485 Quality System - Medical Devices

PROMEPLA
Le Copori
9 avenue ALbert II
MONACO
98000
MONACO

Order Information

Account Number:	104501
Synergy Health Sales Part Reference:	1127465
Customer Reference Number:	426/PLANER - 16/06/2021
Product Description:	496618P - PRODUITS CULTURE CELLULAIRE
Validation Reference:	S12029202
Quantity Received:	3
Customer Minimum Specification kGy:	25.0
Customer Maximum Specification kGy:	45.0
Customer Unit Lot/Batch Number:	G005723

Irradiation Data

Date and Time of Irradiation:	24-JUN-2021 04:09
Calculated Minimum Dose kGy:	29.6
Calculated Maximum Dose kGy:	41.5

Irradiation Release Authorised By Synergy Health Marseille SAS, a STERIS Company

Processing Site: M.I.N. 712 - Arnavaux, , Marseille Cedex 14, 13323 Phone No: +33 (0) 4 91 214 214

Registered Office: M.I.N. 712 - Arnavaux, 13323 Marseille Cedex 14, FRANCE

N° TVA: FR59 343 092 540

Détail de la commande: 426

Référence	Ref client	Quantité	Nb carton	Num lot	N° palette
CLIENT: PLANER PLC					
OOG00199SR148		480	10	G005723	R210075
OOG00199SR148		480	10	G005723	R210076
OOG00199SR148		480	10	G005723	R210077
		1440	30		



<http://www.steris-ast.com>

Certificate of Irradiation

Date Issued: 10-Aug-2021

FR01S12593320-1-1

This is to certify that Synergy Health Marseille, a STERIS Company has where appropriate delivered an irradiation process in accordance with the current certified standards:

EN ISO 11137-1 Sterilisation of Health Care Products
EN ISO 13485 Quality System - Medical Devices

PROMEPLA
Le Copori
9 avenue ALbert II
MONACO
98000
MONACO

Order Information

Account Number:	104501
Synergy Health Sales Part Reference:	1127465
Customer Reference Number:	435/PLANER - 29/07/2021
Product Description:	496618P - PRODUITS CULTURE CELLULAIRE
Validation Reference:	S12029202
Quantity Received:	2
Customer Minimum Specification kGy:	25.0
Customer Maximum Specification kGy:	45.0
Customer Unit Lot/Batch Number:	G005725

Irradiation Data

Date and Time of Irradlation:	10-AUG-2021 05:13
Calculated Minimum Dose kGy:	29.6
Calculated Maximum Dose kGy:	41.9

Irradiation Release Authorised By Synergy Health Marseille SAS, a STERIS Company

Processing Site: M.I.N. 712 - Arnavaux, , Marseille Cedex 14, 13323 Phone No: +33 (0) 4 91 214 214

Registered Office: M.I.N. 712 - Arnavaux, 13323 Marseille Cedex 14, FRANCE

N° TVA: FR59 343 092 540

Détail de la commande: 435

Référence	Ref client	Quantité	Nb carton	Num lot	N° palette
CLIENT: PLANER PLC					
OOG00199SR148		480	10	G005725	R210105
OOG00199SR148		480	10	G005725	R210106
		960	20		



<http://www.steris-ast.com>

Certificate of Irradiation

Date Issued: 26-Jul-2021

FR01S12585763-1-1

This is to certify that Synergy Health Marseille, a STERIS Company has where appropriate delivered an irradiation process in accordance with the current certified standards:

EN ISO 11137-1 Sterilisation of Health Care Products
EN ISO 13485 Quality System - Medical Devices

PROMEPLA
Le Copori
9 avenue ALbert II
MONACO
98000
MONACO

Order Information

Account Number:	104501
Synergy Health Sales Part Reference:	1127465
Customer Reference Number:	428/PLANER - 05/07/2021
Product Description:	496618P - PRODUITS CULTURE CELLULAIRE
Validation Reference:	S12029202
Quantity Received:	7
Customer Minimum Specification kGy:	25.0
Customer Maximum Specification kGy:	45.0
Customer Unit Lot/Batch Number:	G005723 G005725

Irradiation Data

Date and Time of Irradiation:	25-JUL-2021 12:48
Calculated Minimum Dose kGy:	28.2
Calculated Maximum Dose kGy:	40.3

Irradiation Release Authorised By Synergy Health Marseille SAS, a STERIS Company

Processing Site: M.I.N. 712 - Arnavaux, Marseille Cedex 14, 13323 Phone No: +33 (0) 4 91 214 214

Registered Office: M.I.N. 712 - Arnavaux, 13323 Marseille Cedex 14, FRANCE

N° TVA: FR59 343 092 540

Détail de la commande: 428

Référence	Ref client	Quantité	Nb carton	Num lot	N° palette
CLIENT: PLANER PLC					
OOG00199SR148		480	10	G005723	R210079
OOG00199SR148		480	10	G005723	R210080
OOG00199SR148		480	10	G005725	R210081
OOG00199SR148		480	10	G005725	R210082
OOG00199SR148		480	10	G005725	R210083
OOG00199SR148		480	10	G005725	R210084
OOG00199SR148		480	10	G005725	R210085
		3360	70		

PROMEPLA S.A.M.

9, Avenue Albert II
98000 Monaco

TEL : 00 377 97 98 42 32
FAX : 00 377 92 05 61 50
E-mail : sales@promepia.com
Site : www.promepia.com

Monaco le 18/08/21

Packing List N° BL211857

PLANER LTD
110 WINDMILL ROAD
MIDDLESEX
TW16 7HD SUNBURY-ON-THAMES
ANGLETERRE (UK)

Delivery address

PLANER PLC
110 WINDMILL ROAD
UNIT F, DOLPHIN IND ESTATE
TW16 7HD SUNBURY-ON-THAMES
ANGLETERRE (UK)

Your Ref : PNO46271 MAY Code : CPLANER

Shipment : MONACO LOGISTIQUE
N° Tracking : 80CMX1M20X1M94

Number of Cases : 50
Handling Unit : 5
Total Weight : 500

Part N°	Y/Part N°	Description	Batch N° Compl.	Batch N°	Qty	Still to deliver
---------	-----------	-------------	-----------------	----------	-----	------------------

COG00199SR148	CN101558-1	ASSEMBLED BOTTLE AND CAP - DISPOSABLE HUMIDIFICATION ASSY CUSTOMS 9018390000/ EORI NUMBER: FRO0 00 0905900010/ FRANCE ORIGIN		G005723	2400	
---------------	------------	---	--	---------	------	--

cert: FR201512585763-1
-> 260 pcs
cert: FR201512579241-1-1
-> 1060 pcs

PROMEPLA S.A.M.
9 Avenue Albert II
Immeuble Le Copori
98000 MONACO
Phone number : (+377) 97-98-42 32
Mail : contact@promepia.com

CERTIFICATE OF COMPLIANCE

Bearing all circumstances beyond our control, we hereby certify that all the above-mentioned products bearing a batch number are in conformity with the technical specifications and have not undergone any modifications whatsoever with respect to raw materials, processing or technical characteristics and are in accordance with our NF EN ISO13485 : 2016 Quality System.

EXTRACT OF GENERAL TERMS OF SALES
All our goods travel at the address's own risks.
S.A. Monégasque au capital de 586.420 € - N.I.E : FR29000009059 - S.S.E.E. 4669 B0 1526 - R.C.I 69 S 1239 -

PROMEPLA S.A.M.

9, Avenue Albert II

98000 Monaco

TEL : 00 377 97 98 42 32

FAX : 00 377 92 05 61 50

E-mail : sales@promepia.com

Site : www.promepia.com

Monaco le 18/08/21

Packing List N° BL211858

PLANER LTD
110 WINDMILL ROAD
MIDDLESEX
TW16 7HD SUNBURY-ON-THAMES
ANGLETERRE (UK)

Delivery address

PLANER PLC
110 WINDMILL ROAD
UNIT F, DOLPHIN IND ESTATE
TW16 7HD SUNBURY-ON-THAMES
ANGLETERRE (UK)

Your Ref : PND46271 JULY

Code : CPLANER

Shipment : MONACO LOGISTIQUE
N° Tracking : 80cmx1m20x1m94

Number of Cases : 70
Handling Unit : 7
Total Weight : 700

Part N°

Y/Part N°

Description

Batch N° Compl

Batch N°

Qty

Still to deliver

COGCO199SR148

CN101566-1

ASSEMBLED BOTTLE AND CAP - DISPOSABLE HUMIDIFICATION ASSY
CUSTOMS 9018390000/ EORI NUMBER: FR00 00 0905900010/ FRANCE ORIGIN

G005725

3360

*cert: FR015125852963-1-1
-> 2400000000 lot 5225
cert: FR01512593330-1-1
-> 9600000000*

PROMEPLA S.A.M.

9, Avenue Albert II

Immeuble Le Coponi

98000 MONACO

Phone number : (+377) 97 98 42 32

Mail : contact@promepia.com

CERTIFICATE OF COMPLIANCE

Bearing all circumstances beyond our control, we have by certifying that all the above-mentioned products bearing a batch number are in conformity with the technical specifications and have not undergone any modifications unless ever with respect to raw materials, processing or technical characteristics and are in accordance with our NF EN ISO13485 : 2016 Quality System.

EXTRACT OF GENERAL TERMS OF SALES

All our goods travel at the address's own risks.

PLANNER

PLANNER Ltd, 110 Windmill Road, Sunbury-on-Thames
Middlesex, TW16 7HD, United Kingdom

Tel : 44(0)1932 755000

E-mail : sales@planner.com

Web : www.planner.com

DESCRIPTION : STERILIZED BOTTLE ASSEMBLY

REF CN101568-1 **RX ONLY**

CE 1639

LOT G005725

CATALOGUE NUMBER

BATCH CODE

STERILE R



STERILIZED USING IRRADIATION DO NOT REUSE

DO NOT USE IF
THE PACKAGE
IS DAMAGED



MANUFACTURE



DATE OF
MANUFACTURE
2021-06-28



STORE IN A DARK PLACE AT ROOM TEMPERATURE

USE BY



2023-11-28

Do not use if package is damaged.

Неповищуйте, якщо балоні пошкоджені.

Ma ikke anvendes, hvis emballagen er beskadiget.

Неодговаривает условиям до упаковки повреждена.

Et ikke anvendes, hvis pakningen er beskadiget.

Ne pas utiliser si l'emballage est endommagé.

Inhalt bei Beschädigung der Verpackung nicht verwenden.

Μην το χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά.

non utilizzare se la confezione è danneggiata.

Nie stosować, jeżeli opakowanie zostało uszkodzone.

Não utilizar se a embalagem estiver danificada.

No utilizar si el envase está dañado.

Får inte användas om förpackningen är skadad.

Неповищувати, якщо пошкоджені умовики.

EC REP

Advena Ltd, Tower Business centre, 2nd Fl.,
Tower Street, Swatar, BKR 4013 Malta



PLANER

PLANER Ltd, 110 Windmill Road, Sunbury-on-Thames
Middlesex, TW16 7HD, United Kingdom

Tel : +44(0)1932 755000

E-mail : sales@planer.com

Web : www.planer.com

DESCRIPTION : STERILIZED BOTTLE ASSEMBLY

REF CN101568-1

RX ONLY

CE 1639

LOT G005723

CATALOGUE NUMBER

BATCH CODE

STERILE R



STERILIZED USING IRRADIATION DO NOT REUSE

DO NOT USE IF
THE PACKAGE
IS DAMAGED



MANUFACTURE

DATE OF

MM 2021-05-28



STORE IN A DARK PLACE AT ROOM TEMPERATURE

USE BY
2023-10-28

Do not use if package is damaged.

Nepoužívejte, je-li balení poškozeno.

Må ikke anvendes, hvis emballagen er beskadiget.

Niet gebruiken indien de verpakking beschadigd is.

Ei saa käyttää, jos pakkaus on vahingoittunut.

Ne pas utiliser si l'emballage est endommagé.

entart bei Beschädigung Verpackung nicht verwenden.

Μην το χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιό.

Non utilizzare se la confezione è danneggiata.

Nie stosować, jeżeli opakowanie zostało uszkodzone.

Não utilizar se a embalagem estiver danificada.

No utilizar si el envase está dañado.

Får inte användas om förpackningen är skadad.

Nepoužívat, ak je balenie poškodené.

EC REP

Advena Ltd, Tower Business Centre, 2nd Flr.,
Tower Street, Swalea, BKR 4013 Malta



GRN: GR2183

Supplier	PA169 Promepla Sa	Receipt date	27/08/2021	Created From	Purchase #PN238	Order
Supplier advice note	B:211857 / BL211857	Department				

PO line#	Count items ?	Location	Part#	Store	Bin	Serial/Lot#	Expiration	Qty
4	No	Production	CN101568-1	QS	QS/BOTTLE	G005725		3,360
1	No	Production	CN101568-1	QS	QS/BOTTLE	G005723		2,400